Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with FotoDent® gingiva.

To avoid damage of material quality does not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the FotoDent® gingiva material.

Possible dangers

Please note the safety instructions for processing FotoDent® gingiva:

Hazard statements: H315 Causes skin irritation. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. H411 Toxic to aquatic life with long lasting effects.

Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear protective gloves/eye protection/face protection. P308+P313 IF exposed or concerned: Get medical advice/attention. P333+P313 If skin irritation or rash occurs: Get medical advice/attention. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: Urethane acrylate in HPMA, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide, Pigments

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® gingiva packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Product description

FotoDent® gingiva is a light-curing resin for the production of gingival masks using 385 nm LED based stereolithography. For obligatory post curing of FotoDent® gingiva built parts, the high performance light curing units FotoDent® Flash and PCU LED are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® gingiva we recommend wearing personal protective equipment f. i. suitable gloves, safety glasses etc.
- FotoDent® gingiva should be well shaken approx. 1 hour before use.
- Pour FotoDent® gingiva carefully into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
Rapid Shape printer: 23 ± 3 °C
Asiga Freeform Pro2 UV: 35 ± 3 °C
- Select the product specific parameter settings for your process.
- When building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing
FotoDent® flash 10 minutes + 1 minute gas pre-flow
PCU LED 8 minutes, 80 % light power
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® gingiva. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content.