To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/or color variations of the FotoDent® guide material.

Possible dangers

Please note the safety instructions for processing FotoDent® guide.

Hazard statements: **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects. (REF D35600 only) **H413** May cause long lasting harmful effects to aquatic life. (REF D35650 only)

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. (REF D35650 only) **P363** Wash contaminated clothing before reuse. (REF D35600 only) **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-hydroxyethyl methacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® guide packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® guide 385 nm or 405 nm are light-curing resins for the production of dental surgical drill guides using 385 nm or 405 nm LED based stereolithography. For the post curing of FotoDent® guide built parts, the light curing units PCU LED N₂ or FotoDent® flash are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® guide we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- Carefully pour FotoDent® guide into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:

RapidShape D and DII series	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C
- Select the product specific parameter settings for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing **385 nm**:

PCU LED N ₂	5 minutes, 33 % light output
FotoDent® flash	6 minutes + 1 minute gas pre-flow
- Post-curing **405 nm**:

FotoDent® flash	10 minutes + 1 minute gas pre-flow
-----------------	------------------------------------
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® guide. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® guide.